

Universal Aluminum Crutches

- Made of lightweight, anodized aluminum.
- Reinforced center tube and footpiece provide added strength.
- Push button adjustment makes changing height easy.
- Height adjustment in 1" increments.
- Handgrip position is readily adjusted with wing nuts.
- Underarm and handgrip pads are comfortable and durable.
- Includes laminated color operating instructions.
- 300 lb maximum weight capacity.
- Designed to meet the requirements of HCPCS code: E0114.
- Limited lifetime warranty.

3610M-8	Adult	8 pr/cs
3610LF-8	Adult, Latex-Free	8 pr/cs
3611M-8	Tall	8 pr/cs
3611LF-8	Tall, Latex-Free	8 pr/cs
3612M-8	Youth	8 pr/cs
3612LF-8	Youth, Latex-Free	8 pr/cs
GF3612-1	Youth	1 pr



3610M-8

ITEM	Armpad to Grip Height Range	Armpad to Floor Height Range	Recommended Patient Height Range
3610M-8, 3610LF-8	15" – 19"	45" – 53"	5' 2" – 5' 10"
3611M-8, 3611LF-8	16" – 20"	53" – 61 ½"	5' 10" – 6' 6"
3612M-8, 3612LF-8	9" – 13"	37" – 46"	4' 6" – 5' 2"



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com