

## Adjustable Aluminum Canes

- Contoured handle is designed to provide a comfortable, secure hand hold.
- Attractive aluminum finish resists marring and scratching.
- Easy push-button height adjustment with locking collet for added security.
- Height adjusts from 30" to 39" in 1" increments.
- Available with a choice of grips, either durable vinyl or soft foam, and a choice of colors shown below.
- Includes laminated color operating instructions.
- Designed to meet the requirements of HCPCS code: E0100.
- 250 lb maximum weight capacity.
- Limited lifetime warranty.

### Vinyl Grip

**6220A**      Aluminum      6/cs


### Soft Grip

**6221A**      Aluminum      6/cs

**6222A**      Bronze      6/cs

**6223A**      Black      6/cs

**6224A**      Multi-pack      6/cs (2 each: Aluminum, Bronze, Black)



Adjustable  
Aluminum Cane



**6221A:** Foam Grip



**6220A:** Vinyl Grip



**Graham-Field Health Products**

2935 Northeast Parkway  
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

**1-800-347-5678**

**Fax: 1-800-726-0601**

**www.grahamfield.com**