

Breezy 100



Breezy 100

Simply modular...



Modularity made simple...

One chair, one order form, multiple applications. The Breezy 100 provides a wide choice of backrests, rear wheels, accessories and seating options conforming perfectly to the individual needs of each user.

Attractive steel and composite design ensures that the Breezy 100 is a stylish, reliable and adaptable wheelchair.

Swing in, swing out detachable legrest, flip-up and detachable armrest and folding backrest options, ensure that the Breezy 100 is ideal for client transfers, accessing confined places and compact transportation.



Breezy 105



Breezy 115



Breezy 121



Breezy 112



Breezy 110




Breezy 141

Configurable...


Clients can be positioned via a range of headrests, anatomic cushions & backrests, elevating legrest and height adjustable armrest.



Anatomic seat and abductor wedge (toileting cut out as option), anatomic back and seat belt.



Ergonomic elevating legrest. Swing in, swing out and detachable.



Choice of long or height adjustable armrest. Flip-back and detachable as standard.



Breezy 100

Simply modular...

Breezy 100



370-520 mm




SW + 190 mm




410 mm




1080 mm



510 mm



450 mm



960 mm

Technical data sheet

Frame colour:	Blue, gold, chrome-plated, summer red	Reversing width:	1070 mm
Upholstery colour:	grey	Folded width:	340 mm
Approved user weight:	125 kg	Folded height:	960 mm
Max. weight:	16 kg	Folded length:	820 mm
Weight of heaviest part:	14.7 Kg		

Measures bases on model Breezy 105

For further information on the full specification, options and accessories, please refer to the order form.

All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com or alternatively is available on request in large text.



Kahlbachring 2-4
69254 Malsch/HD
Tel.: +49 (0) 72 53 / 9 80-0
Fax: +49 (0) 72 53 / 9 80-111

www.SunriseMedical.com